



qbio  
quantitative  
biology

# QBIO MASTER PROGRAM

## quantitative biology in practice

# IMAGING BIOLOGICAL SYSTEMS

## Practical : building a confocal

OCTOBER 21

Emmanuel Margeat ([Margeat@cbs.cnrs.fr](mailto:Margeat@cbs.cnrs.fr))

## OBJECTIVES

1. Basic design of a confocal microscope : optics
2. Selections of appropriate fluorescence filters
3. Fluorescence correlation spectroscopy : some reminders
4. Measurements
5. Option 1 : How to get an image using the confocal ?
6. Option 2 : How to build a multicolour microscope ?



## EQUIPEMENTS

1. Laser : THORLABS LDM635
2. Objective : Zeiss Plan Achromat 100X, NA 1.4
3. Detector : Perkin Elmer SPCM AQR-14
4. Fiber for detection : Thorlabs M18L01
5. Pinole : 100 $\mu$ m diameter

## MEASUREMENTS

1. Determine the concentration of a concentrated ( $>\mu\text{M}$ ) solution of Atto655 without using its fluorescence properties
2. Determine the concentration of a diluted ( $\sim\text{nM}$ ) solution of Atto655 using FCS. You can use the concentrated solution as a standard
3. Determine the size of fluorescent beads in solution using FCS. For reference, the diffusion coefficient of the Atto655 dye is  $300\mu\text{m}^2/\text{s}$

## FOR THE NEXT SESSION

1. Prepare a design of the confocal microscope with all its components properly placed
2. Propose a protocol to measure the concentration of the two solutions
3. Propose a protocol to measure the size of the fluorescent beads