

Laser Safety Training

Caroline CLERTÉ

Oct. 2021

outlook

- Brief introduction to lasers

Definition, Characterization, type.

Definition of MPE-limit and NOHD.

- Laser Classifications

Overview and explanation of what is included in each laser class.

- Main hazards

Eyes injuries

Skin injuries

Associated with combustible materials.

- Safety laser Protections

Global


Individual

- Good practices in laser laboratory



Introduction


- **Definition:** Laser (Light **A**mplification **S**timulated **E**mission **R**adiation)
- **Characteritics:**
 - The energy generated by the laser is in ranges from UV (100-400nm), visible(400-750nm) and IR (750nm-1 mm).



- Monochromatic
 - Directional
 - Coherent
-
- Continuous wave (CW) (>0.25s)
 - Pulsed Laser (<0.25s)
 - Q-switched Laser (pulse 5-25 ns high power) (used in dermato)




Principle



- **Active Medium:** - Gaz (Ar; He-Ne; CO₂)
- Solid (Cristal Nd-YAG; Fiber silice)
- Liquid (dye solution)
- Semi conductors or diodes
- **Excitation or pumping :** - Optic (laser, flash lamp,...)
- Electric (HT current), ...

Types of lasers



Maximum Permissible Exposure (MPE)

- Lasers output characteristic (Wavelength, Power, Size and Duration of exposure) determine how much exposure to laser is Hazardous.
- Maximum permissible exposure** (MPE- limit) = greatest exposure that people can tolerate without injury.
- It is expressed in terms of allowable exposure time (s) for a given irradiance (W/cm^2) at a particular wavelength.
- If the calculated value for a specific laser $> \text{MPE} \Rightarrow$ collective and individual **safety protections** & procedures should be put in place.

Table 3-3. MPE for Selected Lasers and Exposure Times
(Reference: ANSI Z136.1-2000)


Laser Type	Wavelength (μm)	MPE (average power density—watts/ cm^2)			
		Exposure time in seconds			
		0.25 s	10 s	600 s	3×10^4 s
CO_2	10.6	—	$0.1 \text{ W}/\text{cm}^2$	—	$0.1 \text{ W}/\text{cm}^2$
Nd:YAG (cw) ^a	1.33	—	5.1×10^{-3}	—	1.6×10^{-3}
Nd:YAG (cw)	1.064	—	5.1×10^{-3}	—	1.6×10^{-3}
Nd:YAG Q-switched ^b	1.064	—	17×10^{-3}	—	2.3×10^{-6}
GaAs (diode)	0.840	—	1.9×10^{-3}	—	610×10^{-6}
InGdAlP (diode)	0.670	2.5×10^{-3}	—	—	—
HeNe	0.633	2.5×10^{-3}	—	293×10^{-6}	17.6×10^{-6}
Krypton	0.647	2.5×10^{-3}	—	364×10^{-6}	28.5×10^{-6}
	0.568	2.5×10^{-3}	—	31×10^{-6}	18.6×10^{-6}
	0.530	2.5×10^{-3}	—	16.7×10^{-6}	1.0×10^{-6}
Argon	0.514	2.5×10^{-3}	—	16.7×10^{-6}	1.0×10^{-6}

In any case, the workers exposure must be minimal and $< \text{MPE}$.



Nominal Occular Hazard Distance (NOHD)

- The Nominal Occular Hazard Distance (NOHD) is the distance at one should stand away from the laser beam source to be **< MPE-limit**.



- Low power + High diverging beam ☐ Small NOHD
- High power + collimated beam ☐ Huge NOHD





Laser Classification: Class 1

- Completely safe under normal use.
- Output power $< 0.4 \text{ mW}$
- Higher power laser fully enclosed and sealed or interlocked can also belong to Class 1
- Below the MPE limit \Rightarrow No protection required



Class 2 Lasers

- Only visible light lasers (400-700 nm).
- Output power 0.4 – 1 mW
- Blink reflect (150-200 ms) protect against this class
- Intentional viewing can lead to eye damage.
- Below the MPE limit => No protection required



Class 3 Lasers

- 3R : - Output power from **1 to 5 mW** in continuous mode; Safe if handle correctly. Can cause eye injury when the energy is collected and put into the eye with binoculars for example.
 - 3B: - Output power from **5 to 500 mW**. **Hazardous**. Eye protection required and must include safety interlock and key-switch.
- Can be **hazardous** under **direct** and **specular reflection** viewing conditions
 - **Not hazardous** for **diffuse reflections**




Class 4 Lasers

- **Dangerous**
- **> 500 mW.**
- must include safety interlock and key-switch.
- **Hazardous** under **direct, specular** and **diffuse** reflections conditions for the **skin** and the **eyes**.
- These lasers can also ignite **fire**



MaiTai Fs Pulsed IR laser. 3 W.



Most powerful laser 10PW.




Summary Laser Classification

TYPE	INJURY OR RISK	RISK LEVEL
Class 1	Low energy levels, not hazardous to skin or eyes. Safe during normal operation.	Low
Class 1M	Safe during normal operation but may cause eye injury if viewed with an optical instrument.	Low-Medium
Class 2	Visible wavelengths only, natural blink response provides eye safety.	Low-Medium
Class 2M	Visible wavelengths only, blink response provides eye safety for unaided viewing	Medium
Class 3R	Transitional zone between safe and hazardous laser products. Direct viewing of	Medium-High
Class 3B	Direct viewing and specular reflections can cause eye injury. Diffuse reflections are	High
Class 4	Can cause severe skin and eye injury through any direct exposure, specular reflections and/or sometimes from diffuse reflections.	Extreme



Main hazards




- Damages depend of **Wavelength**, **Power**, Exposure **Duration** and **Type** (Direct, Diffuse reflection or long term low level).

- Eye Injuries
- Skin Injuries

- Explosions and Fire Hazards
- Electrical Hazards
- Chemical Hazards





Eye injuries




- Photochemical
- Retinal burn



Eye injuries



Diffuse light




Laser beam




UV light

From 180-315 nm



- Photokeratitis & eye sunburn


From 315-400 nm



- Photochemical cataract
- Clouding of the eye lens.






IR (780-1400 nm)



- IR Cannot be seen but will focus on the back of the retina causing **cataract** and **retinal burn**.




IR (1.4 μm - 1 mm)




Aqueous flare; cataract and corneal burn.


Impact




cataract



Skin injuries



















Skin injuries




- Damage can range from surface or dermis **skin burn, sunburn, increased of skin aging, or increase of risk of skin cancer (UV)**
- **Photosensitive reactions** can occur (VIS range)




Summary

Risk Class	Eye Direct laser beam vision	Eye Direct laser beam vision with optics	Eye diffuse beam vision	Skin	Fire
1					
1M					
2					
2M					
3R					
3B					
4					

 Hazardous if $t > 0.25$ s

 Could be Hazardous

 Hazardous

 Fire ignition
Hazardous




Safety controls

- **Engineering controls** involved design features applied to laser, laser beam or laser environment that restrict exposure. (eg: beam shutter, beam attenuator, protective housing, interlocks)
- **Administrative control measures** involve procedures and information (eg: Laser Safety Training, Posting warning signs and labels, Nomination of a Laser safety officer, publishing operating procedures, entry limitations)
- **Environment protective equipment** : Laser lab should be equipped with controlled access, safety signs, safety protections (cover, curtains,...) and safety notices.
- **Personal protective equipment**: is worn by personnel using or being in vicinity of lasers. This control includes protective eyewear, gloves and clothing.






Environment safety

- **Laser room:** controlled access, Non reflecting floor, Matt & light paint on the walls, 500 Lux lights, protective curtains, protective covers.
- **Control units:** Key command for laser ON, Emergency stop button, Laser safety signs light on when laser is on, Interlocked keys.




Environment safety

- **Security unit:** Emission alarm, protective cover, shutters, laser power meter, class identification, danger symbol




Laser room safely equipped



Laser safety Glasses

Laser Glasses Selection Guide



European Norm: **EN 207** & **CE-mark**




OD Specs (ANSI Z136)

180 to 400 nm, OD = 6+


720 to 1090 nm, OD = 5+

750 to 1064 nm, OD = 7+

LG9 Optical Density



Entering Laser Room Good practices



- ✓ Remove **all reflective objects** and jewelry
- ✓ If the **Laser sign** is **ON**, **knock on the door** and wait for a response
- ✓ Wear appropriate individual protections (**safety glasses, labcoat**)
- ✓ **Do not turn on the light** if you don't know if an experiment is ongoing.







Buildup & Alignment Good practices



- ✓ Use appropriate laser **safety glasses**.
- ✓ Work with good **light** in room to minimize pupil size.
- ✓ Work with **lowest laser power** as possible.
- ✓ **NEVER** look at the laser **directly**.
- ✓ Use **card** or **camera** to Follow and align the laser beam.



Buildup & Alignment Good practices



- ✓ **Control the propagation** of laser beam (Laser Guarding)
- ✓ Keep the laser beam path in **horizontal plan below eye level**.
- ✓ Always **Fix optic components on optical table**.
- ✓ **Avoid all specular and diffuse reflections.**



Lasers sources @ CBS29

ROOM	Wavelengths	Power mW	Class	Which setup
1A14	LSM800 : 405/488/561/640	CW UV-visible 100	1	Airy scan
	Lumencor Celesta	1000 mW	4	Airy scan
1A14	CW - UV 405 nm	150	3b	MFM Antoine
	CW- visible 488 nm	50	3b	MFM Antoine
	CW- visible 561 nm	150	3b	MFM Antoine
	CW - visible 642 nm	110	3b	MFM Antoine
	CW - IR 785 nm	100	4	MFM Antoine
1A10	CW - IR 850 nm	CW IR 100 mW at 850 nm	Class 4	AFM Super resolution PEM
				AFM Super resolution PEM
	CW-UV-Visible 405;488;568;634 nm	CW 100 mW	Class 3b	AFM Super resolution PEM
				AFM Super resolution PEM
	CW - IR 635 nm	CW laser diode 2 mW	Class 3b	AFM Super resolution PEM
1A11	CW laser diode 6 mW 660nm (51nanoFI-S Shafter + Kirchhoff)	CW laser diode 6 mW	Class 3b	AFM Super resolution PEM
	CW laser diode 30 mW 630 nm (48FI-5-630 Shafter + Kirchhoff)	CW laser diode 30 mW	Class 3b	AFM Super resolution PEM

- 8 lasers class 3B
- 3 lasers class 4



Lasers sources @ CBS60

ROOM	Wavelengths	Power mW	Class	Which setup
OB06	CW - UV 405 nm	100	3b	RAMM
	CW- visible 488	100	3b	RAMM
	CW- visible 561	50	3b	RAMM
	CW- visible 642	100	3b	RAMM
	CW - IR 785	100	4	RAMM
OB08	CW - UV 405	100	3b	PALM
	CW- visible 488	100	3b	PALM
	CW- visible 561	100	3b	PALM
	CW- visible 642	100	3b	PALM
	CW - IR 785	100	4	PALM
	Pulsed IR - 700-1020 nm	3 W à 800nm	4	2PE FCS
	Pulsed and CW 405 nm	50	3b	STED-FCS
	Pulsed and CW 485 nm	50	3b	STED-FCS
	Pulsed and CW 595 nm	5	3b	STED-FCS
	Pulsed and CW 640 nm	50	3b	STED-FCS
OB10	CW- visible 450	>100	3b	MFM+Single Mol Manu
	CW- visible 488	>100	3b	MFM+Single Mol Manu
	CW- visible 514	>100	3b	MFM+Single Mol Manu
	CW- visible 532	500	4	MFM+Single Mol Manu
	CW- visible 640	100	3b	MFM+Single Mol Manu
	CW - IR 785	100	4	MFM+Single Mol Manu
	CW- visible 633	17	3b	2 spots FCS Manouk
OB13	CW- visible 660	100	3b	AFM-Lucas
OB02	Pulsed IR 690-1300 nm	1,46 W à 900 nm	4	2-photon IPAM
OB02	Pulsed IR 680-1080 nm	3 W à 800 nm	4	2-photon IPAM

- 18 lasers class 3B
- 7 lasers class 4




Main causes of Laser injuries

- ❑ **80%** of Eye injuries are due to **specular reflections**.
- ❑ **Inadequate Eye protection**: 60% due to Laser safety glasses not worn.
- ❑ 80% accidents take place during **alignment**.
- ❑ **Bad Control** the **propagation** of laser beam.
- ❑ 75% injured person have **>5 years expertise**.
- ❑ **Bad training, no information or sensibilization**.
- ❑ **No respect of safety procedures**.



Victims Laser Accidents



Laser accident Statistic-professional activity



Actions when accident occur

- 1) Injured person should be take care of **immediately**.
 - **Lay down the victim on the floor.**
- 2) Contact **Hospital ophtalmologist emergency**. **04 67 33 77 90**
Address: Hospital Gui de Chauliac, 80 avenue Augustin Fliche Montpellier.
 - Ophtalmologic control the very day and 10-15 days after.
- 3) **Report** any accident on the register book.



There is no second chance!!

DOCTOR FUN



Peer pressure in the laser lab





Optional slides





EXPOSITION DIRECTE

FAISCEAU DIRECT





FAISCEAU FOCALISÉ

Lentille biconvexe = cas de l'objectif




FAISCEAU DIVERGENT

Lentille biconcave




EXPOSITION À UNE RÉFLEXION SPÉCULAIRE

FAISCEAU DIRECT



Surface polie et plane


FAISCEAU DIVERGENT



Faisceau divergent

Surface polie et convexe

FAISCEAU FOCALISÉ




Faisceau refocalisé


Surface polie et concave

EXPOSITION À UNE RÉFLEXION DIFFUSE

DIFFUSION UNIFORME (Diffuseur parfait)



DIFFUSION NON UNIFORME (En pratique)



Optical Density and Transmission

$$E_x = E_0 10^{-OD} \quad \longrightarrow \quad T = 10^{-OD}$$

- **Transmission** is defined as the ration of E_x / E_0
- If OD of the laser safety eyewear is = 3, the transmission will be of 10^{-3} or 0.1%

Table 3-6. Transmission at Various OD Values

OD	Transmission		%Transmission
0	1	1.0	100%
1	10^{-1}	0.1	10%
2	10^{-2}	0.01	1%
3	10^{-3}	0.001	0.1%
4	10^{-4}	0.0001	0.01%
5	10^{-5}	0.00001	0.001%
6	10^{-6}	0.000001	0.0001%
7	10^{-7}	0.0000001	0.00001%

